

UC60-5-GD-2

Medical Storage Cabinet




Cart shown with optional trays, baskets, catheter slide hooks, and electronic keyless locks

Dimensions (WxHxD): 42.322" x 70.25" x 28.532"


Dimensions include bumpers

41.322"

28.532"



70.25"



Product Description

UNICAS UC60-5-GD-2 double-column medical storage cabinet with tempered glass double doors is designed for efficient organization and quick access to medical supplies.

Features

5" medical grade casters

- (2) Top brake casters (front)
- (2) Swivel casters (back)

Tempered glass doors

Magnetic door latch mechanism

- Conveniently keep cabinet doors securely closed

Optional electronic keyless locks

Customizable storage

- Optional trays, baskets, and accessories to meet your storage needs

Polypropylene bumpers

Clear anodized aluminum uprights and doors

Robust zinc plated steel base frame

Plastic ABS interior panels

- Ability to angle storage trays and baskets for ease of access