

# Medical Procedure Cart

TKWT5100-DB-NH



Model Shown: TKWT5100-DB-NH

## Product Description

The UNICAS medical procedure cart is built around our customizable basket and tray system which is engineered to streamline the efficiency of healthcare procedures. This versatile cart is equipped with precision and convenience in mind, offering a robust solution for medical professionals.

### Features

#### Organized Storage

- The cart offers built-in divisible drawers for optimal organization of supplies

#### Mobility

- With smooth-rolling 5" casters this cart ensures effortless mobility

#### Durable Construction

- Crafted from high-quality materials, our Medical Procedure Cart is built to withstand the demands of a medical environment

#### Customization Options

- Tailor the cart to meet specific needs with customizable drawer configurations

#### Security

- Wing, key, and keyless electronic lock options are available to provide secure storage for sensitive items

#### Enhanced Stability

- Our outrigger stabilizing frame provides enhanced stability, significantly reducing the risk of tipping over

#### 5" medical grade casters

- (2) Top brake casters (front)
- (2) Swivel casters (back)

**Dimensions (WxHxD):** 32" x 45" x 25"

